

seca 959

Digital chair scale
with swivelling arm
and foot rests III



- Patients easy to move from bed to chair due to swivelling arm and foot rests.
- Smooth running wheels.
- BMI function.
- Clear control element.
- Version with PC interface.
- Including carrying case for mains adapter.

seca 959:

Consideration, comfort and safety.

How can one weigh bedridden patients or people with an impaired sense of balance as gently as possible? The answer is a seca chair scale as this can be positioned directly next to the patient, saving the patient – and the carer – unnecessary effort.



Swivelling arm and foot rests for easy transfer to scale.

Consideration and comfort – these are particularly prominent aspects of this stable model with its circular steel tube frame. The armrests which can be folded back make the unimpeded transfer of the patient from chair or bed to the chair scale possible. The foot rests can be folded upwards and pushed under the seat. The lockable rear wheels ensure that it is safe to sit down and stand up.

Clearly laid out control element.

Positioned on the back of the chair scale seca 959 and thus easy to operate is the display and control element with its large LCD display. All clearly laid out functions are easy to control with four large keys. Practical: the carrying case seca 471 provided for storing the mains adapter directly on the scale.

Version with PC interface.

One version is fitted with an interface on the bottom side of the display element so that the results can be transferred directly to a PC. This interface is especially helpful when using the seca 959 in the dialysis sector. The carrying case seca 471 for the mains adapter is included in delivery.

Effortless mobility.

The seca 959 has four smooth running pivoting wheels making use possible even in the narrowest of spaces.

Precise results for all cases.

As the zero point of the weighing function can change due to external influences – especially in the case of such a highly sensitive scale, it is automatically checked and aligned electronically once every second. And to ensure that the seca 959 determines the precise weight of restless patients, it is fitted with an effective damping system. An additional function: the BMI value and therefore the nutritional condition of the person weighed can be ascertained automatically at the push of a button.



The practical positioning of the display allows easy operation.



To make sitting down and standing up easier, the arm and footrests can be folded away.

seca 959

Technical Data

- Capacity: 200 kg
- Graduation: 50 g < 150 kg > 100 g
- Dimensions (WxHxD):
563 x 906 x 970 mm
- Weight: 24.5 kg
- Power supply:
Mains adapter / rechargeable batteries
- Functions: TARE, Hold, BMI,
automatic switch-off
- Approval class:

